

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

Manufacturer and Model	Sonntek P/N	Compatible with OEM P/N	Lamp Type	List Price
ABI (Applied Biosystems)/Kratos/Persp.Bios./Life Tech.				
120A,130A	LH-98	2450-0128	D2	\$450
270A-HT	LH-98		D2	\$450
Analyst	LH-13		D2	\$570
Analyst	LH-13X2		D2	\$750
Analyst	L-850		W	\$205
Biocad 700E,20,60,250,Sprint	LH-40	P5-1085-00-5020	D2	\$555
Biocad 20, 60, 250	LH-40	P5-1085-00-5020	D2	\$555
Integral 250,700E	LH-40	PS-1085-00-5020	D2	\$555
Peptide Synthesizer 433A	LH-98	N2920149	D2	\$450
Prosice 492,494,HT,CLC,C	LH-98	2900-0484	D2	\$450
AGILENT				
1037,1047A	L-610	01037-60002	W	\$385
1046	L-890	2140-0549	P/XE	DISCONTINUED
SPECTROPHOTOMETERS				
8450A	LH-16	08450-60106	D2	\$555
8451A	LH-11	08451-60100	D2	\$555
8452A	LH-12	08452-60104	D2	\$555
8453A,E, 8454	LH-15	2140-0605	D2	\$595
8453A, 8454	LH-15X2	2140-0605	D2	\$780
8453A,E, 8454	L-850	G1103-60001	W	\$205
1040 SERIES				
1040A,1040M,79881A	LH-10	79880-60002	D2	\$520
1050 SERIES				
79853A	L-295	79853-60002	D2	\$420
79853C,79854A,G1306A	LH-10	79883-60002	D2	\$520
1090 SERIES				
79880A,79881A	LH-10	79883-60002	D2	\$520
SERIES I,II	LH-10	79883-60002	D2	\$520
1100/1120/1200/1220 SERIES				
G1314A,B,C,D,E,F; G7114 A,B	LH-14	G1314-60100	D2	\$500
G1315A,B,C,D,E G1365A,B,C,D	LH-13	2140-0590	D2	\$570
G1315A,B,C,D,E G1365A,B,C,D	LH-13X2	2140-0590	D2	\$750
G1315A,B,C,D,E G1365A,B,C,D	L-850	G1103-60001	W	\$205
G1321A	L-870	2140-0600	P/XE	\$785
HP3D/CE	LH-14	2140-0585	D2	\$500
G1314 D,E,F	LH-14	G1314-60101	D2	\$500
G1315C,D, G1365C,D	LH-13X2	2140-0820	D2	\$750
1260,1290 INFINITY SERIES				
G4212 A/B DAD, G7117A/B/C/D, (8 PIN)	LH-19X2	5190-0917	D2	\$750
G4212 A/B DAD, G7117A/B/C/D, (8 PIN)	LH-19	5190-0917	D2	\$570
G1315A,B,C MWD/DAD, G7115A, G7165A	LH-13	G2140-0590	D2	\$570
G1315A,B,C MWD/DAD, G7115A, G7165A	LH-13X2	G2140-0590	D2	\$750
G1315A,B,C MWD/DAD, G7115A, G7165A	L-850	G1103-60001	W	\$205
CAPILLARY ELECTROPHORESIS				
G1600CE,1601CE,G1600AX	LH-14	2140-0585	D2	\$500
G7100 (8-PIN)	LH-19X2	5190-0917	D2	\$750
G4212 A/B DAD (8 PIN)	LH-19X2	5190-0917	D2	\$750
AKTA (SEE GE HEALTHCARE)				
ALLTECH (SEE GRACE)				
ALTEX (SEE BECKMAN)				
AMERSHAM (SEE GE HEALTHCARE)				
AMINCO (SEE THERMO SCIENTIFIC)				
ANALYTIK JENA				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

M500	LH-95	820-60021-0	D2	INQUIRE
S5 (UV),S10 (UV/VIS),S100, S300, S600	LH-95	820-60021-0	D2	INQUIRE
SPEKOL 1200	LH-95	820-60021-0	D2	INQUIRE
SPEKORD 40, 50	LH-95	820-60021-0	D2	INQUIRE
SPEKORD 200, 205, 210, 250	LH-95	820-60021-0	D2	INQUIRE
APPLIED BIOSYSTEMS (SEE ABI)				
APPLIED SCIENCE/ROCHE (SEE ROCHE/APPLIED SCIENCE)				
ATI (SEE THERMO SCIENTIFIC)				
ARGONAUT TECHNOLOGIES (SEE BIOTAGE)				
BAS				
FL45	LH-910	MF4526	XE	\$630
UV116A,UV200	LH-40	MF4650	D2	\$555
BARNSTEAD/TURNER/METERTECH (SEE THERMO SCIENTIFIC)				
BAUSCH & LOMB (SEE THERMO SCIENTIFIC)				
BECKMAN/ALTEX/DANAHER				
20,24,26,34,35,36	L-250	585699	D2	\$390
20,24,25,26,3435,36	L-860	564333	W	\$45
25	L-200		D2	\$395
163	L-205	947029	D2	\$400
166	LHA-70	538706	D2 Lamp Ass'y	\$670
166	LH-70		D2 lamp Only	\$555
168	LHA-75	538711	D2 Lamp Ass'y	\$670
168	LH-75		D2 lamp Only	\$555
5200,ACTA	L-250	585699	D2	\$390
5240	L-840	873947	W	\$40
AD340	L-275		W	\$45
DB,DBG,DBGT,DK,DU	L-305	92680	D2	\$395
DB,DBG,DK-2A,DU2	L-305	92680	D2	\$395
DU-6,7,8,20,30,37,40,50	L-250	585699	D2	\$390
DU-6,7,8	L-840	873947	W	\$40
DU-20,30,40,50	L-835	514259	W	\$30
DU-60,62,63,64,65,68	LH-71	596791	D2	\$670
DU-60,62,63,64,65,68	L-835	514259	W	\$30
DU-70	LH-73	523253	D2	\$670
DU-500,520,530	LH-110	517660	D2	\$555
DU-500,520,530	L-730	517659	W	DISCONTINUED
DU-600,640,640B,650,800	LH-72	514366	D2	\$635
DU-600,640,640B,650,800	L-835	514259	W	\$30
DU-700,720,730	LH-78	A23792	D2	\$650
DU-700,720,730		A23778	W	
DU-7000 SERIES	LH-72	514366	D2	\$635
DU-7000 SERIES	L-835	514259	W	\$30
P/ACE 2000,2050,2100,2200,2210	LH-77	359058	D2	\$650
BECKMAN (continued)				
P/ACE 5000/477170,5010/477171	LH-76	142205	D2	\$680
P/ACE 5500/477172,5510/477173	LH-76	142205	D2	\$680
P/ACE MDQ	LH-76	144667	D2	\$680
P/ACE 800	LH-76	144667	D2	\$680
PROTEOME LAB PF20-166	LHA-70	538706	D2 Lamp Ass'y	\$670
PROTEOME LAB PF20-166	LH-70		D2 Lamp Only	\$555
BEIJING PURKINJE GENERAL INSTRUMENT CO/GENTECH				
BEIJING RALEIGH ANALYTICAL INSTR. CORP (BRAIC)				
UV 1100,9100,1200,1300,1400	L-265	TU-1901	D2	\$385
VIS-7220,9100,1200,1300,1400	L-165		W	\$40
BIOCHROM/HARVARD BIOSCIENCE/GE/PHARMACIA				\$0

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

LIBRA S12	* LH-975		D2/W	\$675
LIBRA S12	L-935	80-2106-16	W	\$30
LIBRA S21,22	L-990		P/XE	\$775
ULTROSPEC 1000 E, I	* LH-975	80-2104-56	D2/W	\$675
ULTROSPEC 2000,3000,4000		80-2104-56		INQUIRE
BIO-RAD				
BIOLOGIC DUO FLOW/QUADTEC				
750-0200	* L-25	750-0216	HG	\$555
750-1300, 760-1300, 760-1330	LH-25		D2	\$555
750-1300, 760-1300	A4073		W Lamp Ass'y	\$260
1706,1806A	L-295	167-0159	D2	\$410
1790, 1801, BIO-DIMENSION	LH-40	167-0905	D2	\$555
BIOFOCUS,MICROSAMPLER 100	LH-41	148-6008	D2	\$555
REMEDI	LH-115	195-7602	D2	\$650
SMARTSPEC 3000	L-60	170-2501	P/XE	\$670
788-4009				INQUIRE
ICYCLER IQ,IQ5	L-24	170-8756	W	INQUIRE
MYIQ,MYIQ2	L-24	170-8756	W	INQUIRE
BIOTAGE/ARGONAUT				
HORIZON	L-735	09427	HG	\$540
ISOLERA SYSTEMS				
ISO-1SF,EF,4SF,EF	L-735	09427	HG	\$540
ISO-1SV,SW,EV,EW,4SV,SW,EW	LH-60	09830	D2	\$670
SP1,SP4,SPX	LH-60	09830	D2	\$670
KILOPREP(100G)	LH-40		D2	\$555
FLASHMASTER II	LH-22		D2	\$555
MXF SERIES				INQUIRE
BIO-TEK INSTRUMENTS, INC. /AGILENT				
EL SERIES				INQUIRE
ELX-808	L-1055		W	\$105
ELX-800, ELX-800WAM	L-1060		W	\$65
ELX800NB, ELX800UV	L-1065		W	INQUIRE
FLX-800T,TB,TBI,TBID,TBX	L-1035		W	INQUIRE
SYNERGY 2/FI/FP	L-1035		W	INQUIRE
SYNERGY H4, HI/FI	L-1035		W	INQUIRE
SYNERGY SERIES			P/XE	INQUIRE
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				
BMG LABTECH				
INQUIRE				
BUCHI				
028062/028063, 028067/028068 C-635 C-630	LH-22	37401	D2	\$555
C-640	LH-60	1105813	D2	\$670
C-640		11058132	W	INQUIRE
PREPCHROM, C-650, C-700	LH-60	11058132	D2	\$670
				\$0
BUCK (SEE CECIL/JENWAY)				
BLC-20	LH-94	7200001	D2	\$555
CAMAG				
TLC II SCANNER	L-285		D2	\$405
CARLO ERBA/THERMO SCIENTIFIC				
	L-280	24439A	D2	\$370
				\$0
CARY/VARIAN				
1,3,4,5E,100,300,400,500	L-230	56-1012-15	D2	\$455
1,3,4,5E,100,300	L-930	56-10013900	W	\$35

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

50	L-1045		P/XE	<i>DISCONTINUED</i>
400	L-110	56-10013900	W	\$40
				\$0
219	L-205	56-101215-00	D2	\$395
219	L-925		W	\$40
500	L-1015	56-10013900	W	\$40
4000, 5000, 6000	* LH-240	110713990	D2	\$500
CECIL (BUCK)				
CE 1020,2020,2040,2041,2501,2502	L-115		D2	\$400
CE 1020,2020,2040,2041,2501,2502	L-845		W	\$55
CE 1200, 1220	LH-20		D2	\$580
				\$0
CE 3021, CE 3031	L-115		D2	\$400
CE 4000 thru 9000 SERIES	L-115		D2	\$400
CE 4000 thru 9000 SERIES	L-845		W	\$55
				\$0
COLE-PALMER				
01953-00	LH-40	H-01953-51	D2	\$555
COMETRO TECHNOLOGY LTD				
6000 PUW				INQUIRE
DIONEX/LC PACKING/GYNKOTEK/THERMO				
<i>(SEE THERMO SCIENTIFIC)</i>				
D-STAR				
DDW-10,DLC-20,DVW-10	LH-94	025-0078	D2	\$555
				\$0
DUPONT				
850	L-280		D2	\$380
8800	L-280		D2	\$380
				\$0
DYAX (SEE BIOTAGE)				
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

ECOM				
LCD2073	LH-60	22590001	D2	\$670
SAPPHIRE,OPAL,TOPAZ	LH-60	22590001	D2	\$670
FLASH 06, FLASH 10	LH-60	22590001	D2	\$670
LCD2071.2, LCD2071.3	L-735	70150001	HG	\$540
EG & G (SEE PERKIN-ELMER)				
EMD MILLIPORE/MERCK MILLIPORE (SEE HITACHI)				
L2455	LH-56	8902430	D2	\$450
NOVA 30A,60A				INQUIRE
PHARO 100,300				INQUIRE
MILLIPORE WATER SYSTEM			HG	\$425
ERMA/ERC				
7511, 7512, 7515, 7522, 7530, 7535	L-610	8811550	W	\$355
ESA				
520, 522	LH-40	70-1924	D2	\$555
FINSTRUMENTS				
INQUIRE				
FISHER (SEE THERMO SCIENTIFIC)				
GBC				
AA LAMPS	L-200	97-0301-00	D2	\$395
CINTRA 906,908,916	L-200	97-0301-00	D2	\$395
CINTRA 916	L-840		W	\$40
CINTRA 1155, 1255	LH-910		P/XE	\$630
CINTRA 1200	L-280	240-0126-00	D2	\$370
CINTRA 1205,1215,1210	LH-40	1003656	D2	\$555
CINTRA 5100	LH-82		D2	\$555
GE HEALTHCARE/CYTIVA/AMERSHAM BIOSCIENCES				
PHARMACIA/LKB/BIOCHROM				
213,821,402,141				
2151	LH-20		D2	\$580
4050, 4051, 4053, 4050	L-210	80-2084-33	D2	\$410
4050, 4060	L-835	80-2022-94	W	\$30
AKTA SERIES/FPLC/PROTEIN PURIFICATION SYSTEM				
AKTA BASIC, ETAN, EXPLORER 10, 100, 105	LH-990	18-1029-35	P/XE	\$1,905
AKTA FPLC, EXPRESS,CROSS FLOW, PRIME PLUS	L-1010	28-9354-93	HG	\$790
AKTA MICROSYSTEM, PILOT	LH-990	18-1029-35	P/XE	\$1,905
AKTA PURIFIER 10, 100	LH-990	18-1029-35	P/XE	\$1,905
AKTA PURIFIER 10, 100	L-1010	28-9354-93	HG	\$790
AKTA PURE 25, 150	LH-880	**NEW**	P/XE	\$1,440
AVANT 25, 150	LH-880	**NEW**	P/XE	\$1,440
GENEQUANT II	* LH-975	80-2104-56	D2/W	\$675
GENEQUANT 100, 1300				
GENEQUANT PRO	* LH-975	80-2104-56	D2/W	\$675
NOVASPEC II	L-835		W	\$30
NOVA SPEC III, PLUS			W	
RNA/DNA CALCULATOR	* LH-975	80-2115-33	D2/W	\$675
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

GE HEALTHCARE/CYTIVA (continued)				
OLIGO PILOT	LH-990		P/XE	\$1,905
NANOVUE SERIES				
			P/XE	
ULTROSPEC SERIES				
ULTROSPEC II,III,PLUS	L-210	80-2084-33	D2	\$410
ULTROSPEC II,III,PLUS	L-835	80-2022-94	W	\$30
ULTROSPEC 500 PRO	* LH-125	80-2106-17	D2/W	\$765
ULTROSPEC 500 PRO	L-935	80-2106-16	W	\$30
ULTROSPEC 1000 E, 1000 I	* LH-975	80-2109-11	D2/W	\$675
ULTROSPEC 1000 E	L-935	80-2106-16	W	\$30
ULTROSPEC 1100 PRO	* LH-975	80-2109-11	D2/W	\$675
ULTROSPEC 1100 PRO	L-935	80-2106-16	W	\$30
ULTROSPEC 2000, 3000, 3000 PRO	* LH-125	80-2106-17	D2/W	\$765
ULTROSPEC 2000, 3000	L-935	80-2106-16	W	\$30
ULTROSPEC 2100 PRO, 3100 PRO	L-640	80-2111-29	P/XE	\$565
ULTROSPEC 3300, 4000, 4300	* LH-125	80-2106-17	D2/W	\$765
ULTROSPEC 5300 PRO, 6300 PRO	* LH-125	80-2106-17	D2/W	\$765
UV 900, UVIS 920	LH-990	18-1029-35	P/XE	\$1,905
UV-1	* L-1030	19-3807-01	HG	
	L-1010	28-4042-25	HG	\$790
UV-MII	L-1070	18-0630-01	HG	\$790
UPC 900	L-1010	28-4042-25	HG	\$790
GENTECH (SEE BEIJING)				
GILFORD				
200, 220, 222, 230, 240, 250, 252, 260	L-205	1282X14	D2	\$395
200, 220, 222, 230, 240, 250, 252, 260	L-875	1282X16	W	\$90
2000, 2400, 2600	L-205	1282X14	D2	\$395
2000, 2400, 2600	L-875	1282X16	W	\$90
FLUORO IV	L-801	1452X14	XE	
RESPONSE I	L-210	21429	D2	\$410
RESPONSE I	L-900	1445X42	W	\$100
STASAR I, II, III, SYVA S-3	L-705	1064X19	W	
GILSON				
115, 116, 117	LH-150	100123	D2	\$455
118UV, 119	LH-151	100326	D2	\$455
122	L-800	0902-1511	XE	\$480
151, 152, 155, 156	LH-151	100326	D2	\$455
160, 170, 171, 172				
GRACE/ALLTECH/BUCHI				
200, 201, 205	LH-40	2037077	D2	\$555
200	LH-760	2037079	W	\$350
305	LH-910	304144	P/XE	\$630
REVELERIS UV1, UV2	LH-40	2037077	D2	\$555
REVELERIS X2	LH-22	2119201	D2	\$555
GROTON				
FD100			P/XE	DISCONTINUED
GYNKOTEK (SEE DIONEX)				
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

HACH				
DR/2700,2800,3800,3900			W	INQUIRE
DR/2000, 2010	L-745	4664700	W	\$145
DR/4000U,DR/4000V	L-730	4951600	W	DISCONTINUED
DR/4000U	LH-110	4951300	D2	\$555
48000, 48100	L-730	4951600	W	DISCONTINUED
DR/5000, 6000	LH-78	A23792	D2	\$650
DR/5000, 6000		A23778	W	
2100N, 2100AN LAB TURBIDIMETERS	L-785		W	DISCONTINUED
57806-18			W	
HARVARD BIOSCIENCE (SEE BIOCHROM)				
HITACHI				
100-10,20,30,40,60,80	L-280		D2	\$370
101, 102, 111	L-280		D2	\$370
124	L-210		D2	\$410
150-20	LH-51	ANO-0424-AS	D2	\$615
200, 300, 330, 340	LH-51	ANO-0424-AS	D2	\$615
5440	L-800	650-1500	XE	\$480
655A-22, -23	LH-51	ANO-0424-AS	D2	\$615
F-1000,1050,1080,1200,2000,2500	L-800	650-1500	XE	\$480
F-3000,3010,4000,4010,4500,7000	L-800	650-1500	XE	\$480
GENE SPEC I	LH-58		D2	
L-2400,2420,2450,2455	LH-56	890-2430	D2	\$450
L-2420,2450,2455	L-835	890-2527	W	\$30
L-2480, 2485	L-800	650-1500	XE	\$480
L-2500, 3000, 3100, 3200	LH-51	ANO-0424-AS	D2	\$615
L-4000, 4500	LH-56	890-2430	D2	\$450
L-4200,4250	LH-51	ANO-0424-AS	D2	\$615
L-4200,4250	L-815	885-1200-AS	W	\$100
L-5440	L-800	650-1500	XE	\$480
L-7000,7300,7400,7401,7420,7450A,7455	LH-56	890-2430	D2	\$450
L-7480,7485,7490	L-800	650-1500	XE	\$480
U-0080D	LH-56	890-2430	D2	\$450
U-1100, 1500	LH-51	ANO-0424-AS	D2	\$615
U-1100	L-815	885-1200-AS	W	\$100
U-1800,1900	LH-57	2J11500	D2	\$615
U-1800	L-935	122-2301	W	\$30
U-1900	L-815	885-1200-AS	W	\$100
U-2000,2001,2010	LH-51	ANO-0424-AS	D2	\$615
U-2000,2001,2010	L-815	885-1200-AS	W	\$100
U-2800,2810,2900,2910	LH-57	2J11500	D2	\$615
U-2800, 2810	L-935	122-2301	W	\$30
U-2900, 2910	L-815	885-1200-AS	W	\$100
U-3000,3010,3110,3210	LH-51	ANO-0424-AS	D2	\$615
U-3000,3010,3110	L-815	885-1200-AS	W	\$100
U-3210	LH-55	DISCONTINUED	D2	
U-3300,3310,3400,3410	LH-51	ANO-0424-AS	D2	\$615

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

HITACHI(continued)				
U-3300,3310	L-815	885-1200-AS	W	\$100
U-3500	LH-56	890-2430	D2	\$450
U-3501, 4001, 4100	LH-51	ANO-0424-AS	D2	\$615
U-4100	L-815	885-1200-AS	W	\$100
HNU				
HORIBA JOBIN YVON (SEE JY)				
APSA-370, APOA-370			P/XE	
IBM				
9420				
9523, 9563, 9565	L-205	4780947	D2	\$395
IMPLEN				
P30,31,33,34,36			P/XE	
INTERCHIM				
450	LH-60	G09630	D2	\$670
PURIFLASH MODELS	LH-60	G09630	D2	\$670
IONICS BUSINESS GROUP (SEE SIEVERS)				
ISCO/TELEDYNE				
205	LH-40	250-1400-00	D2	\$555
3850, V4/1840	L-205	60-3144-002	D2	\$395
COMBIFLASH COMPANION,4X,XL,TS	LH-61	605248001	D2	\$555
COMBIFLASH RF+ PURLON	LH-62	605238123	D2	\$555
COMBIFLASH RF 200(RF200I), RF75	LH-62	605238123	D2	\$555
COMBIFLASH SYSTEM SQ16X	L-205	68-1844-142	D2	\$395
COMBIFLASH RF 150	LH-62	605238123	D2	\$555
COMBIFLASH RF & LUMEN	LH-62	605238123	D2	\$555
JASCO				
870A, -875A	L-295	3063	D2	\$410
7850	L-265	0035	D2	\$385
CD-995, -2095	L790		XE-HG	<i>DISCONTINUED</i>
FP-750,-777,-820,-821,-920,-921	L-800	3022	XE	\$480
FP-1520,-1570,-1575,-2020,-3120	L-800	3023	XE	\$480
FP-6200,-6300,-6500,-6600,-8000	L-800	3022	XE	\$480
J-810 CD, -815 CD	L-960	3027	XE	\$605
J-720 (150W)	L-960		XE	\$605
J-720 (450W)	L-1000		XE	
MD 2010, MD 4010	L-265	6678-H-103A		\$385
OR-2090	L-790		XE-HG	<i>DISCONTINUED</i>
UV-970A, -975A	L-295	3066	D2	\$410
UV-970B,-975B,-975C		3067	D2	
UV-1570,-1575,-2070,-2075,-2077		5330.0091	D2	
UV-3070, -3075			D2	
UVIDEC 100,100-II,-III,-IV,-V,-VI	L-200		D2	\$395
UVIDEC 320, 340, 420, 430	L-200		D2	\$395
UVIDEC 505, 510, 610, 650, 710	L-200		D2	\$395
V-500, -530, -550, -560	L-265	0035	D2	\$385
V-530, -550, -560, -570	L-835	0020	W	\$30
V-630 BIO	L-265	0035	D2	\$385
V-650, 650 PC & IRM TYPE, -660	L-265	0020	D2	\$385

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>JASCO (continued)</u>				
V-650	L-835	0020	W	\$30
V-7100	L-265		D2	\$385
X-LC MODELS				
JENWAY/BIBBY SCIENTIFIC				
6105		615-005	D2	
6105		012.075	W	
6300, 6320D, 6400	L-835	012.075	W	\$30
6405, 6505	LH-650		D2	DISCONTINUED
6405, 6500, 6505	L-835	012.075	W	\$30
6700, 7300, 7305			W	
7310, 7315			P/XE	
AQUANOVA, GENOVA NANO/PLUS		012.094	P/XE	
JY (JOBIN YVON)/SPEX/HORIBA				
1934D		96166	P/XE	
FLUOROLOG-2/ F111AI, F112AI		1906	XE	
FLUOROLOG-2/F111XI, F112XI	L-1000	1907	XE	
FLUOROLOG-2/F1T11I, F212I, F222I	L-1000	1907	XE	
FLUOROG-3, FLUOROLOG TAU-3	L-855		XE	DISCONTINUED
FLUOROMAX-2,-3,-4	L-955		XE	DISCONTINUED
SYSTEM 5000	L-60		P/XE	\$670
KNAUER/SONNTEK				
<u>AZURA</u>				
2.1S, 2.1S FIBER OPTIC	A5193	A5193	D2	\$925
2.1 L UVD, MWD, DAD	A5193	A5193	D2	\$925
6.1 L	AZL01	AZL01	D2	
6.1 L	AZL02	AZL02	W	
<u>BLUE SHADOW</u>				
10D, 20D, 40D, 50D	A5193	A5193	D2	\$925
<u>SMARTLINE</u>				
200	A4142	A4142	HG	
2500, 2600	A4071	A4071	D2	\$810
2500, 2600	LH-25	<i>A4071 equivalent</i>	D2	\$555
2500	A4073	A4073	W	
2600	A4072	A4072	W	
2520	A5193	A5193	D2	\$925
2550	A5194	A5194	D2	\$1,290
2550	A5195	A5195	W	\$440
2800, 2850	A4447V1	A4447V1	D2	\$1,290
2800, 2850	A4448	A4448	W	\$440
<u>WELLCHROM</u>				
200	A4142	A4142	HG	
2000, 2001	A4071	A4071	D2	\$810
2000, 2001	LH-25	<i>A4071 equivalent</i>	D2	\$555
2001	A4073	A4073	W	
2500, 2501, 2600, 2650, 2800	A4071	A4071	D2	\$810
2500, 2501, 2600, 2650, 2800 (A4071)	LH-25	<i>A4071 equivalent</i>	D2	\$555
2500, 2501, 2600	A4073	A4073	W	
2800	A4448	A4448	W	\$440

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>KNAUER/SONNTEK (CONTINUED)</u>				
<u>NON-SERIES</u>				
A59200, A59201	A59210	A59210	XE	\$2,730
A59200, A59201	L-795	A59210 equivalent	XE	
6700, 8700, 9700, D871	A0877		D2	DISCONTINUED
6700, 8700, 9700, D871	LH-20		D2	\$580
KONTRON/BIOTEK/NORTH STAR SCIENTIFIC				
440	LHA-34		Lamp Ass'y	\$820
440	LH-34		D2 Lamp only	\$540
720	L-205		D2	\$395
810, 820, 860, 900, 930, 940	L-205		D2	\$395
820, 860, 900, 930, 940	L-110		W	\$40
UVIKON	L-205		D2	\$395
ALL SPECTROPHOTOMETERS	L-205		D2	\$395
KRATOS/SCHOEFFEL (ALSO SEE ABI /PERKIN-ELMER)				
FS950		2450-0099	HG	
FS970, 980	LH-98	SFA-101	D2	\$450
SF757, SF759A, SF769, SF770, SF773	LH-98	N2920149	D2	\$450
SF783, SF783A, 785, 785A	LH-98	N2920149	D2	\$450
VISION	LH-40	5-2466-18	D2	\$555
SF769, SF770, FS970	LHA-1	SFA-101	Lamp Ass'y	\$815
ALL OTHERS	LHA-3		Lamp Ass'y	\$565
LABALLIANCE/SCIENTIFIC SYSTEMS (SSI)				
200, 201, 500, 504, 506, V520, 525	LH-40	201305	D2	\$555
200, 201, 500, 504, 506, V520, 525	LH-760	201305	W	\$350
<u>1200/1500 SERIES</u>				
AGG1314B	LH-14		D2	\$500
AGG1315D	LH-13		D2	\$570
AGG1321A	L-870		P/XE	\$785
LC PACKINGS (SEE DIONEX)				
LDC (SEE THERMO SCIENTIFIC)				
LIFE (SEE ABI)				
LINEAR/OEM/PRIVATE LABEL				
<i>(SEE ALLTECH, BAS, BIO-RAD, COLE-PALMER, DIONEX, FISHER, GBC, JSCO, LABALLIANCE</i>				
<i>LDC/THERMO, MICHROM BIORESOURCES, PERSEPTIVE BIOSYSTEMS, RAININ</i>				
<i>SPECTRA-PHYSICS, THERMO, THERMO SCIENTIFIC, VARIAN)</i>				
LKB (SEE GE HEALTHCARE)				
MCPHERSON				
FL-750/A, FL-750/B	LH-101	FL-750-003	XE	
FL-750/A, FL-750/B	LH-100	FL-750-005	D2	
FL-750/A, FL-750/B	LH-102	FL-750-017	XE-HG	
MERCK (SEE HITACHI)				
METERTECH (SEE BARNSTEAD/TURNER)				
METROHM				
MICHROM BIORESOURCES/BRUKER				
UV-200	LH-40	200-37000-00	D2	\$555
MD1, 2, 3, 4	LH-25	220/25000/00	D2	\$555
MICROTECH SCIENTIFIC, INC.				
LINEAR MT200/201, 205, 200CE	LH-40		D2	\$555

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

MILTON-ROY (SEE THERMO SCIENTIFIC/BAUSCH&LOMB)				
MOLECULAR DEVICES				
ANALYST AD, HT		26-002	XE	
EMAX, EMAX PLUS			W	
FLEX STATION, I, II, II 384	L-870	2100-0014	P/XE	\$785
FLEX STATION III	L-775	4200-0025	P/XE	\$920
OPTIMAX, SPECTRAMAX PLUS	L-60	4200-0044	P/XE	\$670
SPECTRAMAX 190, 340 PC, 384 PLUS	L-60	4200-0044	P/XE	\$670
SPECTRAMAX 250	L-870	2100-0014	P/XE	\$785
SPECTRAMAX 340	L-65	52100-0008	P/XE	
SPECTRAMAX M2/M2E,M3,M4,M5/M5E	L-775	4200-0025	P/XE	\$920
SPECTRAMAX ANALYST GT	L-870	2100-0014	P/XE	\$785
SPECTRAMAX GEMINI,EM,XPS,XS	L-870	2100-0014	P/XE	\$785
SPECTRAMAX i3, PARADIGM				
SPECTRAMAX L				
THERMO MAX	L-120	0310-1001R	W	\$105
UV MAX, VMAX			W	
NICOLET (SEE THERMO SCIENTIFIC/ATI)				
OCEAN OPTICS				
D1000		AISD UV-1		
DH2000		DH2000-BD		
PX-2		PX-2B	P/XE	
OLIS				
RSM-1000	L-1000		X/E	
PARAYTEC				
CLARITY 17, 620, 1000				
TDA200 (100D)	L-870			\$785
PERKIN-ELMER/EG&G				
HPLC UV DETECTORS				
785A	LH-98	N2920149	D2	\$450
INTEGRAL 4000	L-210	N2351125	D2	\$410
LC-55, 55B, 65, 75, 85B, 95	L-210	N2351125	D2	\$410
LC-135, 235	L-210	N2351285	D2	\$410
LC-200	LH-98	N2920149	D2	\$450
LC-200		24500213	W	
LC-90, 250, 290, 295				
LC-480	LH-92	B0160917	D2	\$520
FLEXAR				
UV/VIS LC	LH-98		D2	\$450
UHPLC,PDA, FLUORESCENCE				
200/DAD, PHASE I (SHORTER VERSION FOR OLDER MODELS)	LH-97		D2	
200/DAD PHASE II, III	LH-99		D2	\$570
200/DAD	L-725		W	DISCONTINUED
200/FLUORO	L-800		XE	\$480
200/UV-VIS/785A	LH-98		D2	\$450
200/UV-VIS/785A			W	
PDA PLUS DETECTOR	L-211	09290900	D2	\$650
FLUOROMETERS				
200, 203, 204, 512, 650 SERIES	L-800	03501646	XE	\$480

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>PERKIN-ELMER (continued)</u>				
HTS 7000			P/XE	
LS-3,4,5,30,40,50,50B,55		04969486	P/XE	
LC-240			P/XE	
MPF-2A, 3, 4	L-800		XE	\$480
MPF-43, 43A, 44, 44A, 44B	L-800		XE	\$480
<u>SPECTROPHOTOMETERS</u>				
46, 55, 57	L-200	C0550505	D2	\$395
55			W	
220,320,330,340	LH-51	02390354	D2	\$615
3840	L-200	C6760011	D2	\$395
550, 551, 555, 559, 559A	LH-92	102014	D2	\$520
552A,553,559A,570,571,572,575,576	L-200	C0550505	D2	\$395
556, 557	LH-51	02390354	D2	\$615
BIO	LH-92	B0160917	D2	\$520
BIO	LH-845	B0114620	W	\$125
EZ150, EZ201, EZ210	LH-51	C6957529	D2	\$615
EZ150, EZ201, EZ210	L-840	C6957528	W	\$40
LAMBDA 1, 3, 3B	L-280	C0550505	D2	\$370
LAMBDA 2	LH-92	B0160917	D2	\$520
LAMBDA 1, 2, 3B	LH-845	B0114620	W	\$125
LAMBDA 4, 6	LH-93	C6880055	D2	\$500
LAMBDA 4, 6	LH-845	B0114620	W	\$125
LAMBDA 5,7,9,10,11,12,14,15	LH-92	B0160917	D2	\$520
LAMBDA 5,7,9,10,11,12,14,15	LH-845	B0114620	W	\$125
LAMBDA 16,17,18,19,20,25,35,40,45	LH-92	B0160917	D2	\$520
LAMBDA 16,17,18,19,20,25,35,40,45	LH-845	B0114620	W	\$125
LAMBDA 650,750,800,850,900,950,1050	LH-92	B0160917	D2	\$520
LAMBDA 650,750,800,850,900,950,1050	LH-845	B0114620	W	\$125
LAMBDA BIO, BIO PLUS,XLS, XLS PLUS			P/XE	
<u>ATOMIC ABSORPTION SPECTROMETERS</u>				
100, 300 ANALYSIS	L-335	N0370119	D2	
600, 700, 800 ANALYSIS	L-405	B0148615	D2	\$395
3100, 3110, 3300 ANALYSIS	L-335	N0370119	D2	
4000 ANALYSIS	L-310	00470385	D2	\$395
272, 280, 360, 370, 372, 373, 380, 460, 560	L-310	00570194	D2	\$395
1100, 1100B, 2100	L-405	B0148615	D2	\$395
1100, 1100B, 2100	L-405	B0148615	D2	\$395
2280, 2380, 3030B	L-310	00570194	D2	\$395
4100	L-405	B0148615	D2	\$395
5000, 5100	L-310	00470385	D2	\$395
<u>PLATE READERS</u>				
2100,2103,2104 ENVISION			P/XE	
<u>POLARIMETERS</u>				
240, 241, 243, 341, 343				
PERSEPTIVE BIOSYSTEMS/PERKIN-ELMER/EG&G/APPLIED BIOSYSTEMS (SEE ABI)				
PHARMACIA (SEE GE HEALTHCARE/BIOCHROM)				
PHOTON TECHNOLOGY INT'L (PTI)				
QM-2, 4, 6, 7, 8, 14 /2005	L-970		XE	DISCONTINUED
PICODROP, LTD.				
PICO 100, 200			P/XE	
POLYMER LABS/VARIAN/AGILENT				
PL-ELS 1000	L-840		W	\$40
POSTNOVA (SEE SYCAM)				
S3210, S3240				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

3211, 3212, 3241				
PYE UNICAM (SEE THERMO SCIENTIFIC/UNICAM)				
RAININ (SEE VARIAN)				
ROCHE/APPLIED SCIENCE				
LIGHT CYCLER LC 480	L-1040	04686136001	XE/R	INQUIRE
LIGHT CYCLER LC 96	L-1040	04686136001	XE/R	
ROSEMOUNT ANALYTICAL/EMERSON PROCESS MGT				
890	L-330	654843	D2	\$395
SCHOEFFEL (SEE KRATOS/ABI/APPLIED BIOSYSTEMS)				
SCIENTIFIC COMPUTING INT'L				
FILMTEK				
4000, 4000 EM-DUV		10295		
SCIENTIFIC SYSTEMS (SSI) (SEE LABALLIANCE)				
SCINCO				
S-3100, S-4100				
SEDEX				
55, 65, 75, 75C	L-935	45800	W	\$30
55LT, 60LT, 75LT, 754				
80LT, 85LT			LED	
SEQUOIA TURNER				
340, 360, 690			W	
SHIMADZU				
<u>HPLC UV DETECTORS</u>				
LC2010AHT, LC2010CHT	LH-84	228-37401-00	D2	\$555
SP-4	L-305		D2	\$395
SPD-2A,2AV,6A,7A,7AV	L-305	062-65056-03	D2	\$395
SPD 6AV	L-100		W	DISCONTINUED
SPD-10A, 10AV, 10AVP, 10AVVP	LH-82	228-34016-02	D2	\$555
SPD-10AV, 10AVVP	L-830	670-14602-00	W	\$200
SPD-10AVI	LH-82	228-34016-02	D2	\$555
SPD-20A, 20AV	LH-82	228-34016-02	D2	\$555
SPD-20AV	L-830	670-14602-00	W	\$200
SPD-M6A, SPD-M10A		670-14502-00	D2	
SPD-M6A, SPD-M10A	L-935	062-65005-00	W	\$30
SPD-M10AV		670-14604-00	D2	
SPD-M10AV	L-245	228-34410-91	W	\$30
SPD-M10AVP	LH-88	228-34016-00	D2	\$555
SPD-M10AVP		228-33189-91	W	
SPD-M20A, M20AV	LH-87	228-3401600	D2	\$555
SPD-M20A	L-245	228-34410-91	W	\$30
SPD-M20AV	L-830	670-14602-00	W	\$200
SPD-M30A				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

SHIMADZU (continued)				
ELSD-LTII(2), ELSD-LT				
FLUOROMETERS				
RF-10A, 10AXL		228-34216-00	XE	INQUIRE
RF-10AVP				INQUIRE
RF-20A, 20AXS	* L-795	228-51511-95	XE	
RF-1501, 5000, 5301	L-800	200-81500-01	XE	\$480
RF-530, 535, 540, 551	L-800	204-03133-01	XE	\$480
SPECTROPHOTOMETERS				
BIOSPEC-MINI, BIOSPEC-NANO			P/XE	
BIOSPEC 1601, MULTISPEC 1501	L-265	062-65055-05	D2	\$385
BIOSPEC 1601	L-110	062-65004-06	W	\$40
MULTISPEC 1501	L-935	062-65005-00	W	\$30
CS-920, 930, 9301	L-265	062-65055-05	D2	\$385
CS-930		200-39053-00	W	
SERIES 160	L-110	062-65004-06	W	\$40
UV-120, 150, 160A, 160U, 163	L-265	062-65055-05	D2	\$385
UV-120, 160A, 160U	L-110	062-65004-06	W	\$40
UV-120-01, UV-120-02	L-265	062-65055-05	D2	\$385
UV-120-01, UV-120-02			W	
UV-140-02, UV-150-02	L-265	062-65055-05	D2	\$385
UV-240, 250, 260, 265	L-265	062-65055-05	D2	\$385
UV-240, 260, 265	L-110	062-65004-06	W	\$40
UV-1200, 1201, 1201PC, 1203, 1240	L-265	062-65055-05	D2	\$385
UV-1200, 1201, 1203, 1240, 1240V	L-245	062-65005-00	W	\$30
UV-1501, 1600, 1601, 1601PC, 1603, 1650	L-265	062-65055-05	D2	\$385
UV-1501, 1600, 1601, 1601PC, 1603, 1650	L-110	062-65004-06	W	\$40
UV-1700, 1800	L-265	062-65055-05	D2	\$385
UV-1700, 1800	L-245	062-65005-00	W	\$30
UV-2100, 2101, 2401PC, 2450	L-265	062-65055-05	D2	\$385
UV-2101PC, 2401, 2450	L-110	062-65004-06	W	\$40
UV-2501PC, 2550, 2600, 2700	L-265	062-65055-05	D2	\$385
UV-2501PC, 2550, 2600, 2700	L-110	062-65004-06	W	\$40
UV-3000				
UV-3100, 3100PC, 3100PC	L-265	062-65055-05	D2	\$385
UV-3101PC	L-110	062-65004-06	W	\$40
UV-3600, 3700	L-265	062-65055-05	D2	\$385
UV-3600	L-110	062-65004-06	W	\$40
ATOMIC ABSORPTION SPECTROMETERS				
AA-625,630,640,660,670,680U	L-265	062-65055-05	D2	\$385
AA-6200, -6300	L-265	062-65055-05	D2	\$385
AA-6501, 6601, 6650, 6701, 6800				
SHODEX				
RI71,72,73,74	L-610	SD100900	W	\$385
RI101		8814315	W	
SICKMAIHAKE				
GM 31, 32	L-1005		D2	\$370
	* LH-280		D2 Lamp Ass'y	\$580
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

SIEVERS/IONICS INSTR. BUSINESS GROUP/GE ANALYTICAL INST.				
400, 800	L-700	AER00020-01	HG	\$320
900, 5310C, PPT	L-605	CARK35001	HG	\$320
SKALAR				
SA6711	LH-20		D2	\$580
SA6712, 6750	LH-22		D2	\$555
6850	LH-28		D2	
SONNTEK (SEE KNAUER)				
SPECTRA-PHYSICS (SEE THERMO SCIENTIFIC)				
SPECTRONICS (SEE THERMO SCIENTIFIC)				
SPECTRONICS CORP				
TRANSILLUMINATORS				
UV-B		BLE1T158	HG	
SPEX (SEE JY /JOBIN YVON)				
SPECTROVISION (SEE GROTON/DIONEX)				
STRATAGENE/AGILENT				
MX 3000, 3005 P	L-985		W	DISCONTINUED
SUNCHROM				
SPECTRAFLOW 500	LH-40		D2	\$555
SPECTRAFLOW 501	LH-25	969-200.201	D2	\$555
SPECTRAFLOW 600 DAD				
SYKAM/LASERCHROM				
S3210 UV, S3210 DAD		4010002	D2	
S3210 UV, S3210 DAD		4010001	W	
TECAN				
GENIOS PLUS/GENIOUS PRO	L-1025	10122179	P/XE	
INFINITE F200, NANOQUANT	L-1025	10122179	P/XE	
M200	L-640		P/XE	\$565
INFINITE F500, M1000				
SAFIRE II, ULTRA EVOLUTION	L-870			\$785
SPECTRO FLUOR PLUS	L-1025	10122179	P/XE	
SUNRISE		30000854/3709008	W	
ULTRA	L-870		P/XE	\$785
TELEDYNE/TEKMAR INSTRUMENTS				
PHOENIX 8000		14-7023-024	D2	
TOX ROUND GRID LAMP		512092	D2	
OZONE GENERATOR LAMP 7			D2	
THERMO SCIENTIFIC/DIONEX/LC PACKING/GYNKOTEK				
AMINCO				
AB-2 (SERIES 2)			XE	
SLM 2000	L-280		D2	\$370
SLM 4800, 8000, 8100	L-855		XE	DISCONTINUED
SPF-500		FA182	XE	
ATI				
240, 250	LHA-34		D2 Lamp Ass'y	\$820
240, 250	LH-34		Lamp Only	\$540
BARNSTEAD/TURNER/METERTECH				
110, 111, 112	L-185	110-850	HG	\$110
SP-850	L-180		W	\$40

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>THERMO (continued)</u>				
SP-830+ , SP-870+	L-190		W	DISCONTINUED
SP-890, SP 8001, M965				
<u>BAUSCH & LOMB/MILTON ROY</u>				
AQUAMATE	LH-54	9423UV90002E	D2	\$615
AQUAMATE, AQUAMATE UV	L-910	9423UV90004E	W	
BIOMATE 3, 3-S			P/XE	
BIOMATE 5	LH-54	9423UV90002E	D2	\$615
BIOMATE 5	L-910	9423UV90004E	W	
DE 50	L-225	343424	D2	\$370
EPSILON	L-915	335423	W	\$55
EVOLUTION 60				
EVOLUTION 6000	LH-54	9423UV90002E	D2	\$615
EVOLUTION 6000	L-910	9423UV90004E	W	
GENESYS 2, 5	* LH-105	343425	D2	\$555
GENESYS 2, 5	L-915	335423	W	\$55
GENESYS 8	LH-54	9423UV90002E	D2	\$615
GENESYS 8	* LH-105	343425	W	\$555
GENESYS 10, 20, 30	L-915	335423	W	\$55
SPECTRONIC 20, 20 PLUS, 20D, 20D PLUS	L-905	333384	W	\$35
SPECTRONIC 21,21DV,21MV,21UVD,21UVM	L-945	333307-602	W	\$125
SPECTRONIC 21D, 21DUV	L-225	343424	D2	\$370
SPECTRONIC 70, 100	L-920	333323-104	W	\$30
SPECTRONIC 200	L-915	335423	W	\$55
SPECTRONIC 301, 401	L-915	335423	W	\$55
SPECTRONIC 501, 601			W	
SPECTRONIC 601,700,701,710	L-225	343424	D2	\$370
SPECTRONIC 700,701,710	L-920	333323-104	W	\$30
SPECTRONIC 1001,1201,2000	L-225	343424	D2	\$370
SPECTRONIC 1001, 1201	L-810	335024	W	\$150
SPECTRONIC 2000	L-920	333323-104	W	\$30
SPECTRONIC 3000	L-270		D2	DISCONTINUED
SPECTRONIC 3000			W	
<u>DIONEX</u>				
ULTIMATE MWD 3000, 3000 (RS)	LH-45	6074.1110	D2	\$685
ULTIMATE DAD 3000				
ULTIMATE DAD 3000, MWD 3000	LH-236		W	DISCONTINUED
AD25, PDA-100, UltiMate PDA-3000	LH-3	6080.2001	D2	
VWD 3000,VWD3100 (RS),VWD 3400(RS)	LH-45	6074.1110	D2	\$685
DAD-3000(RS),MWD-3000(RS),VWD 3000,VWD3100 (RS),VWD 3400(RS)	LH-236		W	DISCONTINUED
VDM-2	L-235	040651	D2	\$380
ULT-UVD, UVD-3000	LH-25	160063	D2	\$555
ULT-UVD	A4073	ULT-DT-H1	W	
RF 2000	L-800	5057.1000	XE	\$480
ENVIRONMENTAL INSTRUMENTS/SO2 ANALYZER				
43A, 43CT, 43CH, 43i(TLE), 580B	L-40		P/XE	DISCONTINUED
FISHER				
20-0344	LH-40	39720	D2	\$555
* MUST RETURN SPENT LAMP (& HOLDER ASSEMBLY)				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>THERMO (continued)</u>				
GYNKOTEK				
AD20, PDA20	LH-5	45230	D2	\$630
AD25,PDA100,PDA2996	LH-3	939016T	D2	
RF1002/2000	L-800	5057.1	XE	\$480
UDM-340S	LH-91	5053.1204	D2	\$490
UVD160,170S,170U,320,340S,340U	LH-91	5053.1204	D2	\$490
<u>LDC-MILTON ROY</u>				
SC100-234	LH-40	9551-0023	D2	\$555
SC100-234	LH-760	9551-0111	W	\$350
SM I, II, III, SM D	L-205	108035	D2	\$395
SM 3000, 31000, 4000	L-205	108035	D2	\$395
SM 3100X, 4100	LH-42		D2	DISCONTINUED
SM 5000	LH-47		D2	DISCONTINUED
UV100-411, UV150-511	LH-40	9551-0023	D2	\$555
UV100-411, UV150-511	LH-760	9551-0111	W	\$350
<u>LINEAR</u>				
100, 200, 203, 204, 205, 215PHD	LH-40	9551-0023	D2	\$555
100, 200, 203, 204, 205	LH-760	9551-0111	W	\$350
304, 305	LH-910	9551-0144	P/XE	\$630
<u>MICROPLATE INSTRUMENTS</u>				
<u>NANODROP</u>				
2000, 2000C, 8000, 3300			P/XE	
1000	L-2010		P/XE	DISCONTINUED
<u>NICOLET</u>				
EVOLUTION 60		335908-1305	P/XE	
EVOLUTION 100, 500, 600	LH-54	9423UV90002E	D2	\$615
EVOLUTION 100	L-910	9423UV90004E	W	
EVOLUTION 500, 600	L-940	9423UV90001E	W	\$60
EVOLUTION 300 BB	L-645	10096001	P/XE	
<u>SPECTRA-PHYSICS/FINNIGAN</u>				
CE500,1000,2000	LH-40	9551-0023	D2	\$555
CE500,1000,2000	LH-760	9551-0111	W	\$350
FL2000,3000	LH-910	9551-0144	P/XE	\$630
SC100,200	LH-40	9551-0023	D2	\$555
SC100,200	LH-760	9551-0111	W	\$350
SF100,200	LH-40	9551-0023	D2	\$555
SF100,200	LH-760	9551-0111	W	\$350
SF770,970	LHA-1		D2 Lamp Ass'y	\$815
SF770,970	LH-1		Lamp Only	\$455
SP8400,8440,8440XR,8450,8490	L-215	A1480-010	D2	\$365
SP8480XR	LHA-3		D2 Lamp Ass'y	\$565
	LH-1		Lamp Only	\$455
	LH-1		Lamp Only	\$455
UV100,150,1000,2000,3000,3000HR	LH-40	9551-0023	D2	\$555
UV100,150,1000,2000,3000	LH-760	9551-0111	W	\$350
<u>SURVEYOR</u>				
6057-60030, 6057-60040 <i>may need adapter (included)</i>	LH-46	108052	D2	\$540
6057-60030, 6057-60040	L-765	803247	W	\$290

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<i>THERMO (continued)</i>				\$0
FL PLUS DETECTOR	L-800	00950-01-00005	XE	\$480
SRVYR-PDA1P, 5P ACCELA, UV1P <i>may need adapter (included)</i>	LH-46	108052	D2	\$540
SRVYR-PDA1P, 5P ACCELA, UV1P	L-765	803247	W	\$290
UV6000LP, UV 8000 <i>may need adapter (included)</i>	LH-46	108052	D2	\$540
UV6000LP, UV 8000	L-765	803247	W	\$290
<i>UNICAM</i>				
UV 1-100,2-100,2-200,2-300,2-400	LH-54	335824	D2	\$615
UV 1-100	L-910	9423UV90004E	W	
UV 1-100,2-100,2-200,2-300,2-400	L-940	9423UV90001E	W	\$60
UV-3-100,3-200,4-100,4-200,4-500	LH-54	\$335,824	D2	\$615
UV-3-100,3-200,4-100,4-200,4-500	L-940	9423UV90001E	W	\$60
UV 300, 500	LH-54	335824	D2	\$615
UV 300, 500, 550	L-940	9423UV90001E	W	\$60
SP8800	L-255		D2	<i>DISCONTINUED</i>
4225	LH-40	113610	D2	\$555
HELIOS ALPHA,BETA,GAMMA	LH-54	335824	D2	\$615
HELIOS ALPHA,BETA,GAMMA	L-950		W	<i>DISCONTINUED</i>
HELIOS EPSILON	L-915	335423-000	W	\$55
VARIOSKAN 1671624	L-965		P/XE	<i>DISCONTINUED</i>
<i>TRANSGENOMIC</i>				
FDX-99-2100	L-800	070030	XE	\$480
WAVE 3500,3500A,3500HT,4500,4500HT, DNA, OLIGO	LH-56	070024	D2	\$450
WAVE 4500 HT-HS, WAVE HS	L-800	070030	XE	\$480
<i>TURNER/BARNSTEAD (SEE THERMO SCIENTIFIC)</i>				
<i>UNICAM (SEE THERMO SCIENTIFIC)</i>				
<i>UNICO</i>				
1100, 1200, 2100, 2102	L-125		W	<i>DISCONTINUED</i>
<i>VARIAN/CARY/RAININ/AGILENT</i>				
<i>CARY</i>				
1, 3, 4, 5, 14, 15, 16, 17	L-230	56-101215-00	D2	\$455
1, 3, 4, 5, 14, 15, 16, 17	L-930	56-100217-00	W	\$35
50, 60	L-1045	110639690	P/XE	
87/8700	LH-20		D2	\$580
100, 210, 219, 220, 300, 400, 500	L-230	56-101215-00	D2	\$455
100, 210, 219, 220, 300	L-930	56-100217-00	W	\$35
400	L-110	56-100217-00	W	\$40
500	L-1015	56-10021700	W	\$40
2050	L-290		D2	<i>DISCONTINUED</i>
2200, 2300, 2400	L-205	56-180001-00	D2	\$395
2550	L-295	00-97874-00	D2	\$410
<i>ECLIPSE</i>			P/XE	
9070, 9075	LH-910	01-900052-01	P/XE	\$630
AA200,240,275,280,300,400,600,640	L-230	56-101215-00	D2	\$455
AA775,800,875,975,1275,1475,3000	L-230	56-101215-00	D2	\$455
DMS70,80,90,100,200,300	L-230	56-101215-00	D2	\$455
FL-1, 2	L-305	03-916156-91	D2	\$395
LC-5000, 5500	L-305	03-916156-91	D2	\$395

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

<u>Varian/Cary (continued)</u>				
PRO STAR SERIES				
PDA-2	LH-82	03-926074-01	D2	\$555
PRO STAR, PRO STAR 320	LH-950	R007200556	P/XE	\$565
PRO STAR 310, 9050 ** RE-USE HOLDER **			D2	
PRO STAR 325, 335	LH-325	0110715400	D2	\$650
PRO STAR 330	LH-56	03-935705-02	D2	\$450
PRO STAR 340, 345	LH-40	R000088515	D2	\$555
PRO STAR 340, 345	LH-760	R000088519	W	\$350
PRO STAR 360	LH-910	19-00052-01	P/XE	\$630
PRO STAR 363	L-800	03-926131-02	XE	\$480
SPECTRA A10,20,30,40,640,880	L-230	56-101215-00	D2	\$455
STAR 3050	L-305	03-916156-91	D2	\$395
UV-1, 1A	LH-950	R007200556	P/XE	\$565
UV50,100,200	L-305	03-916156-91	D2	\$395
UV-C,DII,M	LH-40	88-515	D2	\$555
VISTA 5000,5500	L-305	03-916156-91	D2	\$395
WATERS				
401, 403, 404				
420	L-770	78405	HG	\$35
440, 441	L-820		HG	DISCONTINUED
450	LHA-1		D2 Lamp Ass'y	\$815
	LH-1		Lamp Only	\$455
470, 474	L-800	47447	XE	\$480
480, 481	L-205	99499	D2	\$395
484	LH-30	80357	D2	\$475
486	LH-31	80678	D2	\$475
490E		46660	P/XE	
990, 991, 994	LHA-34	21516	D2 Lamp Ass'y	\$820
996	LH-33	WAT052586	D2	\$555
2420, 2424 ELSD	L-120	201000159	W	\$105
2475	LH-925	201000131	XE	DISCONTINUED
2487, 2488	LH-36X2	WAS081142	D2	\$500
2996	LH-33	WAT052586		\$555
ACQUITY UPLC				
FLR DETECTOR	LH-930	201000193	HG/XE	DISCONTINUED
PDA 2998, TUV 2489, eLambda	LH-38	201000281	D2	\$730
eLAMBDA (thermally managed), PDA 2998, TUV 2489 WITH CHIP	LH-131	700005269	D2	
LC MODULE 1/34551 SYSTEM	LH-31	80678	D2	\$475
YOUNG LIN				
YL 9120, 9160				
ZEISS				
MCS 501				

SONNTEK RESEARCH LAMP GUIDE 2025

CALL: 201-236-9300 EMAIL: sonntekinc@gmail.com www.sonntekinc.com

We are looking for Distributors!
201-236-9300

Ask about Distributor, Quantity and Educational Discounts

LEGEND					
D2	DEUTERIUM				
HG	MERCURY				
L	LAMP ONLY				
LH	LAMP & HOLDER				
LHA	LAMP HOLDER ASSEMBLY				
NA	SODIUM				
P/XE	PULSED XENON				
W	TUNGSTEN				
XE	XENON				
X2	EXTENDED WARRANTY				
ZN	ZINC				

Lamp types and prices subject to change without prior notice. Sonntek assumes no liability for representations made herein. All information is accurate to the best of our knowledge. Majority of light sources listed are genuine Sonntek parts made to meet or exceed original equipment manufacturer's specifications. Sonntek is not affiliated with any companies mentioned; company or product names may be trademarks of their respective owners. 01/01/2025.

© Sonntek, Inc. 2025, All Rights Reserved

www.sonntekinc.com